

**(19) World Intellectual Property
Organization
International Bureau**



Rec'd PCT/PTO 30 DEC 2004

(43) International Publication Date
8 January 2004 (08.01.2004)

PCT

(10) International Publication Number
WO 2004/002357 A1

- (51) **International Patent Classification⁷:** **A61C 8/00,**
A61J 1/90, 1/14

(21) **International Application Number:**
PCT/SE2003/001021

(22) **International Filing Date:** 18 June 2003 (18.06.2003)

(25) **Filing Language:** English

(26) **Publication Language:** English

(30) **Priority Data:**
0202052-7 1 July 2002 (01.07.2002) SE

(71) **Applicant (for all designated States except US):** **NOBEL BIOCARE AB (publ)** [SE/SE]; Box 5190, S-402 26 Göteborg (SE).

(72) **Inventor; and**

(75) **Inventor/Applicant (for US only):** **HALL, Jan** [SE/SE]; Stabbegatan 2A, S-416 80 Göteborg (SE).

(74) **Agent:** **OLSSON, Gunnar;** Nobel Biocare AB (publ), Box 5190, S-402 26 Göteborg (SE).

(81) **Designated States (national):** AE, AG, AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, BZ, CA, CH, CN, CO, CR, CU, CZ, DE, DK, DM, DZ, EC, EE, ES, FI, GB, GD, GE, GH, GM, HR, HU, ID, IL, IN, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MA, MD, MG, MK, MN, MW, MX, MZ, NO, NZ, OM, PH, PL, PT, RO, RU, SC, SD, SE, SG, SK, SL, TJ, TM, TN, TR, TT, TZ, UA, UG, US, UZ, VC, VN, YU, ZA, ZM, ZW.

(84) **Designated States (regional):** ARIPO patent (GH, GM, KE, LS, MW, MZ, SD, SL, SZ, TZ, UG, ZM, ZW), Eurasian patent (AM, AZ, BY, KG, KZ, MD, RU, TJ, TM), European patent (AT, BE, BG, CH, CY, CZ, DE, DK, EE, ES, FI, FR, GB, GR, HU, IE, IT, LU, MC, NL, PT, RO, SE, SI, SK, TR), OAPI patent (BF, BJ, CF, CG, CI, CM, GA, GN, GQ, GW, ML, MR, NE, SN, TD, TG).

Published:
— with international search report


For two-letter codes and other abbreviations, refer to the "Guidance Notes on Codes and Abbreviations" appearing at the beginning of each regular issue of the PCT Gazette.

Published:

— *with international search report*

For two-letter codes and other abbreviations, refer to the "Guidance Notes on Codes and Abbreviations" appearing at the beginning of each regular issue of the PCT Gazette.

- (54) Title:** ARRANGEMENT AND METHOD FOR PRESERVING THE NEW-BONEFORMING EFFECT OF GROWTH-STIMULATING SUBSTANCES FOR AN IMPLANT PRODUCT



(S7) **Abstract:** In an arrangement and method, use is made of a container for preserving, even for a long time, the new-bone-forming effect of growth-stimulating substances (GS) applied to one or more implant products. The container is arranged so as, dependent on being acted upon, to allow accessibility of the product concerned with applied GS as the time of its use. The container is arranged to enclose the product or the products with applied GS in an environment which is essentially free from air, water and moisture. The container functions as a vacuum container, and the storage times for GS can be extended considerably in relation to previously used arrangements and methods.

WO 2004/002357 A1